


sempermed®

YOUR GLOBAL GLOVE EXPERT™



PRODUCT CATALOG

got gloves?

We Are Your Global GLOVE EXPERT

With MORE THAN 100 YEARS of glove manufacturing experience, we have long been a pioneer within the disposable glove industry. We are able to leverage the latest MANUFACTURING TECHNOLOGY in quality and **innovation** to ensure our gloves achieve optimal **safety**, **reliability**, and **quality**.



GROWING & HARVESTING

True global integration.
From beginning to end, we focus on modern glove manufacturing techniques, research and development of innovative hand protection products and rigorous quality control.



GLOVE EXPERTS

Gloves are our only business. Our product assortment includes latex, nitrile, vinyl, polyethylene, and chloroprene disposable gloves. We have expertise in every phase of the manufacturing process.



WORLD CLASS MANUFACTURING

Setting standards of excellence. Our glove facilities are some of the most sophisticated manufacturing plants in the world. We utilize state-of-the-art technology to achieve the highest levels of quality in our medical and industrial grade gloves.



RESEARCH & DEVELOPMENT

Science never stands still. Neither do we. We are constantly conducting advanced research and development in our ongoing commitment to develop new products that will complement our existing lines.



We Are People Protection.®

Rev08242021

got gloves?

Follow Us!



SCAN ME

sempermed®
YOUR GLOBAL GLOVE EXPERT™



Our Gloves are Disposable. The People they Protect are Not.®
www.sempermedusa.com • Ph 1.800.366.9545

CHOOSING THE RIGHT GLOVE

GLOVE MATERIAL	ATTRIBUTES	ISSUES	PROPERTIES	APPLICATIONS
<p>Chloroprene <i>For maximum protection, change gloves every 30 minutes</i></p> 	<ul style="list-style-type: none"> • Excellent barrier protection • Greater puncture resistance than latex • Comfort: similar to latex • Finish: textured 	<ul style="list-style-type: none"> • Accelerators in glove manufacturing process may cause skin irritation (Less common compared to latex allergies) 	<ul style="list-style-type: none"> • Molecular structure allows wearer to quickly identify a rip, hole or tear • Displays good flame resistance and performs well in a wide range of temperatures 	<p>IDEAL FOR PROLONGED USE & APPROPRIATE FOR MOST APPLICATIONS</p> <p>medical, dental, extreme temperatures, chemical exposure, lab, dermatology, jan/san, safety, and beauty</p>
<p>Nitrile <i>For maximum protection, change gloves every 20 - 30 minutes</i></p> 	<ul style="list-style-type: none"> • Excellent barrier protection • Greater puncture resistance than latex • Comfort: similar to latex • Available in a variety of colors • Finish: textured 	<ul style="list-style-type: none"> • Accelerators in glove manufacturing process may cause skin irritation (Less common compared to latex allergies) 	<ul style="list-style-type: none"> • Molecular structure allows wearer to quickly identify a rip, hole or tear • Less susceptible to material breakdown and leakage than a vinyl or latex glove 	<p>APPROPRIATE FOR MOST APPLICATIONS</p> <p>medical, dental, vet, chemical exposure, lab, jan/san, safety, food services, beauty, and tattoo</p>
<p>Latex <i>For maximum protection, change gloves every 20 - 20 minutes</i></p> 	<ul style="list-style-type: none"> • Excellent barrier performance • Appropriate for most uses • Highly elastic, comfortable fit • Finish: micro-textured and/or fully textured 	<ul style="list-style-type: none"> • Potential for allergies to latex proteins. Reaction can range from sneezing to anaphylactic shock (life-threatening) • Allergies affect 10 to 20% of work-exposed individuals 	<ul style="list-style-type: none"> • Strength and elasticity • Barrier protection against a variety of chemicals • Less susceptible to material breakdown and leakage 	<p>APPROPRIATE FOR MOST APPLICATIONS</p> <p>general patient examinations, dental, wound care, routine oral procedures, lab, food services, beauty, and jan/san</p>
<p>Vinyl <i>For maximum protection, change gloves every 15 minutes or less</i></p> 	<ul style="list-style-type: none"> • Most economical choice • Comfort: not critical • Available in a variety of colors • Finish: smooth 	<ul style="list-style-type: none"> • Cuff has loose fit and is susceptible to material breakdown and leakage • Natural body oils can breakdown the glove material 	<ul style="list-style-type: none"> • Minimal elasticity/flexibility • Typically clear in color, however colored vinyl is widely used in food manufacturing 	<p>INTENDED FOR SHORT-TERM NON-INVASIVE TASKS</p> <p>dispensing medication, non-invasive physical examinations, jan/san, food services, and food manufacturing</p>
<p>Polyethylene <i>For maximum protection, change gloves frequently</i></p> 	<ul style="list-style-type: none"> • Most economical choice for non-medical markets • Finish: smooth 	<ul style="list-style-type: none"> • Baggy, loose fit • Not intended for medical use or stringent chemicals • Very thin material • Intended for very short-term tasks only 	<ul style="list-style-type: none"> • Minimal elasticity/flexibility • Typically clear in color • Low protection • Susceptible to ripping/tearing and leakage 	<p>INTENDED FOR LIGHT DUTY APPLICATIONS</p> <p>jan/san, food services, cleaning</p> <p>* NOT FOR MEDICAL USE</p>



SUSA REORDER #	PRODUCT NAME	GRADE	CHEMO TESTED	ESD TESTED	SKIN FRIENDLY	FENTANYL	CBD OIL	GASTRIC ACID	SPECIAL COATING	GLOVE COLOR
NITRILE										
BKNF	SemperForce® Black	Exam	—	✓	—	✓	✓	—	—	Black
BLNF	SemperForce® Blue	Exam	✓	✓	—	✓	—	—	—	Blue
GSBNF	GripStrong® Black	Industrial	—	—	—	—	—	—	—	Black
GSNF	GripStrong®	Industrial	—	✓	—	—	—	—	—	Black
GSWNF	GripStrong® White	Industrial	—	✓	—	—	—	—	—	White
INIHD	SemperGuard® HD	Industrial	—	✓	—	—	—	—	—	Blue
INIPFT	SemperGuard®	Industrial	—	✓	—	—	—	—	—	Blue
NIPFT	SemperCare®	Exam	✓	✓	—	✓	—	—	—	Blue
ORNF	SemperForce® Orange	Exam	—	—	—	✓	—	—	—	Orange
SMNP	StarMed® PLUS	Exam	✓	—	—	✓	—	—	—	Blue
SMNR	Best Touch® StarMed® Rose	Exam	—	—	—	—	—	—	Aloe	Pink
SMNS	StarMed® Select	Exam	✓	—	—	✓	—	—	—	Blue
SMTN	StarMed® Ultra	Exam	✓	—	—	✓	—	—	—	Blue
SSEC	SemperShield® Extended Cuff	Exam	✓	—	—	✓	—	—	—	Blue
SSSC	SemperShield® Standard Cuff	Exam	—	✓	—	✓	—	—	—	Blue
SUNB	SemperSure® Black	Exam	✓	✓	✓	✓	—	✓	Accelerator-Free	Black
SUNF	SemperSure®	Exam	✓	—	✓	✓	—	✓	Accelerator-Free	Blue
SUNG	SemperSure®	Exam	✓	—	✓	—	—	—	Accelerator-Free	Green
TTNF200	TenderTouch® 200	Exam	✓	—	—	✓	✓	—	—	Blue
TTNF100	TenderTouch® 100	Exam	✓	—	—	✓	✓	—	—	Blue
LATEX										
GSLF	GripStrong®	Industrial	—	—	—	—	—	—	—	Orange
INDPFT	SemperGuard®	Industrial	—	✓	—	—	—	—	—	Orange
INDPS	SemperGuard®	Industrial	—	—	—	—	—	—	Powdered	Orange
SCLT	SemperCare®	Exam	—	—	—	—	—	—	—	Orange
SM	StarMed®	Exam	—	—	—	—	—	—	—	Orange
VINYL										
EVNP	SemperCare® Synthetic	Exam	—	—	—	—	—	—	—	Orange
SCVNP	SemperCare® Stretch	Exam	—	—	—	—	—	—	—	Orange
GSVF	GripStrong®	Industrial	—	—	—	—	—	—	—	Orange
VBLK	SemperGuard® Black	Industrial	—	—	—	—	—	—	—	Black
VP	SemperGuard®	Industrial	—	—	—	—	—	—	Powdered	Orange
VPF	SemperGuard®	Industrial	—	—	—	—	—	—	—	Orange
POLYETHYLENE										
GSPE	GripStrong® Poly	Industrial	—	—	—	—	—	—	—	Orange
HPEF	SemperGuard® FoodSafe Stretch Poly	Industrial	—	—	—	—	—	—	—	Orange

Glove Chemical Resistance Chart*

*This chart is provided as a guideline only. It does not suggest that any Sempermed® glove has been tested for use with the chemical listed. Sempermed® accepts no responsibility for end-user glove selection against a specific risk. Conditions of use and appropriate glove choice are the responsibility of the user.

E = Excellent **G** = Good **F** = Fair **P** = Poor **N** = Not Rated



CHEMICAL	NITRILE (NBR)	NATURAL RUBBER LATEX (NRL)	VINYL (PVC)	CHLOROPRENE (NEOPRENE)
Acetaldehyde	P	G	G	G
Acetic Acid	F	E	E	E
Acetone	P	G	F	P
Alcohol (Butyl)	E	E	E	E
Alcohol (Ethyl)	E	E	G	E
Alcohol (Isopropyl)	E	E	G	E
Ammonium Hydroxide (conc.)	E	G	E	F
Aniline	E	F	F	G
Ammonia	P	G	F	P
Antifreeze (methanol-based)	E	E	P	E
Benzaldehyde	F	F	G	F
Benzene	P	P	F	F
Benzyl Chloride	N	F	P	P
Bleach	P	G	F	P
Brake Cleaner (containing ethanol or hexane)	E	P	P	P
Brake Cleaner, non-chlorinated (containing acetone, n-heptane or xylene)	E	P	P	P
Bromine	N	G	G	G
Calcium Hypochlorite	G	P	G	G
Carbon Disulfide	G	P	F	P
Carbon Tetrachloride	F	P	F	F
CBD Oil	E	N	N	N
Chlorine	E	G	G	G
Chloroacetone	N	F	P	E
Chloroform	P	P	P	F
Chromic Acid	G	P	E	F
Cyclohexane	E	F	P	E
Cytotoxic Chemicals	E	G	P	E
Dibenzyl Ether	N	F	P	G
Dibutyl Phthalate	G	F	P	G
Diesel Fuel	E	P	P	E
Diethanolamine	G	F	E	E
Diethyl Ether	E	F	P	G
Dimethyl Sulfoxide	E	N	N	N
Ethyl Acetate	P	F	F	G
Ethylene Dichloride	P	P	P	F
Ethylene Glycol	E	G	E	G
Fentanyl	E	G	F	E
Fluorine	N	G	G	G
Formaldehyde	E	G	E	E
Formic Acid	E	G	E	E
Gasoline (oil)	E	P	P	E

CHEMICAL	NITRILE (NBR)	NATURAL RUBBER LATEX (NRL)	VINYL (PVC)	CHLOROPRENE (NEOPRENE)
Gastric Acid	G	P	P	G
Grease	G	P	P	G
Glucose	G	P	P	G
Glycerol	E	G	E	G
Hair Dye	G	P	P	N
Hand Sanitizer	E	E	G	E
Hexane	E	P	P	E
Hydrochloric Acid (conc.)	G	P	F	F
Hydrofluoric Acid (30%)	F	P	F	E
Hydrogen Peroxide	G	G	E	G
Insecticides / Pesticides	P	G	F	P
Iodine	G	G	G	E
Kerosene	E	P	P	E
Methylamine	E	G	E	G
Methyl Cellosolve	P	F	P	E
Methyl Chloride	P	P	P	E
Methylene Chloride	P	F	F	F
Methyl Ethyl Ketone	P	F	P	G
Mineral Spirits	P	G	F	P
Monoethanolamine	N	F	E	E
Morpholine	P	F	E	E
Motor Oil	G	P	P	N
Naphthalene	E	G	G	G
Nitric Acid (conc.)	P	P	G	P
Ortho-phthalaldehyde	G	P	P	G
Paint (latex based)	G	P	P	N
Paint (oil-based)	G	P	P	N
Perchloric Acid	F	F	E	G
Peroxide	G	P	P	N
Phenol	F	G	E	G
Phosphoric Acid	E	G	E	E
Potassium Hydroxide (sat.)	E	G	E	G
Propane	E	P	P	N
Sodium Hydroxide	E	G	E	E
Sodium Hypochlorite	E	G	G	E
Sulfuric Acid (conc.)	F	G	G	E
Toluene	P	P	F	F
Trichloroethylene	P	P	F	F
Tricresyl Phosphate	E	P	F	F
Triethanolamine	E	F	E	E
Turpentine	E	P	P	N
Xylene	F	P	N	P

REV04202022



NITRILE Disposable Gloves



SemperGuard[®] HD Nitrile

FEATURES

- Industrial grade heavy duty nitrile
- ESD tested and approved
- Fully textured
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 6-7 mil glove

Order # Size

INIHD102	S
INIHD103	M
INIHD104	L
INIHD105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
72 (12 x 6)



All Sempermed gloves are single use only.



SemperGuard[®] Nitrile

FEATURES

- Industrial grade blue nitrile
- Textured fingertips
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 4-5 mil glove

Order # Size

INIPFT102	S
INIPFT103	M
INIPFT104	L
INIPFT105	XL
INIPFT106	XXL

GLOVE COUNT
100 Gloves/Box (90 XXL)
10 Boxes/Case

PALLET CONFIGURATION
78 (13 x 6)



All Sempermed gloves are single use only.



GRIPSTRONG[®] BLUE NITRILE

FEATURES

- Industrial grade blue nitrile
- Textured fingertips
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove

Order # Size

GSBLN102	S
GSBLN103	M
GSBLN104	L
GSBLN105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
90 (15 x 6)



All Sempermed gloves are single use only.

REV05182023



sempermed[®]
YOUR GLOBAL GLOVE EXPERT[™]

100+ YEARS
EST. 1920



NITRILE Disposable Gloves



GRIPSTRONG® NITRILE

FEATURES

- Industrial grade blue nitrile
- Textured fingertips
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove



Order # Size

GSNF102 S

GSNF103 M

GSNF104 L

GSNF105 XL

GLOVE COUNT

100 Gloves/Box

10 Boxes/Case

PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.



GRIPSTRONG® *WHITE* NITRILE

FEATURES

- Industrial grade white nitrile
- Textured fingertips
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove



Order # Size

GSWNF102 S

GSWNF103 M

GSWNF104 L

GSWNF105 XL

GLOVE COUNT

100 Gloves/Box

10 Boxes/Case

PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.



GRIPSTRONG® BLACK NITRILE

FEATURES

- Industrial grade black nitrile
- Textured fingertips
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove



Order # Size

GSBNF102 S

GSBNF103 M

GSBNF104 L

GSBNF105 XL

GLOVE COUNT

100 Gloves/Box

10 Boxes/Case

PALLET CONFIGURATION

78 (13 x 6)

All Sempermed gloves are single use only.

REV05182023



sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920



NITRILE Disposable Gloves



SEMPERSHIELD®

FEATURES

- Exam grade nitrile - Fully textured
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 6-7 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order # Size

SSSC102	S
SSSC103	M
SSSC104	L
SSSC105	XL
SSSC106	XXL

GLOVE COUNT

100 Gloves/Box (90 XXL)
10 Boxes/Case

PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



SEMPERSHIELD® EXTENDED CUFF

FEATURES

- Exam grade nitrile - Extended cuff
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Tested for use against Fentanyl
- Thickness: Competes with a 5-6 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order # Size

SSEC102	S
SSEC103	M
SSEC104	L
SSEC105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION

78 (13 x 6)

All Sempermed gloves are single use only.



SemperCARE®

FEATURES

- Exam grade nitrile - Fully textured
- ESD tested and approved
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Thickness: Competes with a 4-5 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order # Size

NIPFT101	XS
NIPFT102	S
NIPFT103	M
NIPFT104	L
NIPFT105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.

REV05182023



sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+
YEARS
EST. 1920



NITRILE Disposable Gloves



SEMPERFORCE **BLACK**

FEATURES

- Exam grade black nitrile
- Textured fingertips
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 5-6 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order #	Size
BKNF102	S
BKNF103	M
BKNF104	L
BKNF105	XL
BKNF106	XXL

GLOVE COUNT

100 Gloves/Box (90 XXL)
10 Boxes/Case

PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.



SEMPERFORCE **BLUE**

FEATURES

- Exam grade blue nitrile
- Textured fingertips
- Meets USP 800 requirements for PPE compliance for chemotherapy
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 5-6 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order #	Size
BLNF102	S
BLNF103	M
BLNF104	L
BLNF105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.



SEMPERFORCE **ORANGE**

FEATURES

- Exam grade orange nitrile
- Fully textured
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 5-6 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order #	Size
ORNF102	S
ORNF103	M
ORNF104	L
ORNF105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.

REV05182023



NITRILE Disposable Gloves



StarMed® ULTRA Nitrile Exam

FEATURES

- Exam grade nitrile - Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Thickness: Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order #	Size
SMTN251	XS
SMTN252	S
SMTN253	M
SMTN254	L
SMTN255	XL

GLOVE COUNT

250 Gloves/Box (225 XL)
10 Boxes/Case

PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



StarMed® PLUS Nitrile Exam

FEATURES

- Exam grade nitrile
- Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Thickness: Competes with a 3 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order #	Size
SMNP301	XS
SMNP302	S
SMNP303	M
SMNP304	L
SMNP305	XL

GLOVE COUNT

300 Gloves/Box (250 XL)
10 Boxes/Case

PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



StarMed® Select Nitrile Exam

FEATURES

- Exam grade nitrile
- Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Thickness: Competes with a 3 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order #	Size
SMNS101	XS
SMNS102	S
SMNS103	M
SMNS104	L
SMNS105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case


PALLET CONFIGURATION

64 (16 x 4)

All Sempermed gloves are single use only.

REV05182023





NITRILE Disposable Gloves



Best Touch[®] StarMed[®] ROSE Nitrile Exam with Aloe

FEATURES

- Exam grade nitrile
- Textured fingertips
- Coated with aloe
- Moisturizes dry, chapped hands
- Thickness: Competes with a 3-4 mil glove



Order # Size

SMNR201	XS
SMNR202	S
SMNR203	M
SMNR204	L
SMNR205	XL

GLOVE COUNT

200 Gloves/Box (180 XL)
10 Boxes/Case

PALLET CONFIGURATION

72 (12 x 6)

All Sempermed gloves are single use only.



Tender Touch[®] 200 Nitrile Exam

FEATURES

- Exam grade nitrile - Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order # Size

TTNF201	XS
TTNF202	S
TTNF203	M
TTNF204	L
TTNF205	XL

GLOVE COUNT

200 Gloves/Box (180 XL)
10 Boxes/Case

PALLET CONFIGURATION

72 (12 x 6)

All Sempermed gloves are single use only.



Tender Touch[®] 100 Nitrile Exam

FEATURES

- Exam grade nitrile - Textured fingertips
- Tested for protection against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order # Size

TTNF101	XS
TTNF102	S
TTNF103	M
TTNF104	L
TTNF105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.

REV05182023



Food Handling Compliant



sempermed[®]
YOUR GLOBAL GLOVE EXPERT[™]

100+ YEARS
EST. 1920



NITRILE Disposable Gloves



SemperSure® Black

FEATURES

- Accelerator-Free exam grade nitrile
- Tested for protection against Fentanyl & gastric acid
- Meets USP 800 requirements for PPE compliance for chemotherapy
- ESD tested & approved
- Thickness: Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order # Size

SUNB101	XS
SUNB102	S
SUNB103	M
SUNB104	L
SUNB105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION

64 (16 x 4)

All Sempermed gloves are single use only.



SemperSure®

FEATURES

- Accelerator-Free exam grade nitrile
- Textured fingertips
- Low dermatitis potential
- Tested for protection against Fentanyl & gastric acid
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Thickness: Competes with a 3-4 mil glove

For more information about protection against fentanyl, please contact your sales rep.



Order # Size

SUNF201	XS
SUNF202	S
SUNF203	M
SUNF204	L
SUNF205	XL

GLOVE COUNT

200 Gloves/Box (180 XL)
10 Boxes/Case

PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



SemperSure® Green

FEATURES

- Accelerator-Free exam grade nitrile
- Textured fingertips
- Innovative green technology
- No chlorine used in manufacturing process
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Thickness: Competes with a 3-4 mil glove

Order # Size

SUNG201	XS
SUNG202	S
SUNG203	M
SUNG204	L
SUNG205	XL

GLOVE COUNT

200 Gloves/Box (180 XL)
10 Boxes/Case

PALLET CONFIGURATION

60 (10 x 6)

All Sempermed gloves are single use only.



REV05182023



sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920

5 REASONS to Choose

SemperSure® Green



2 MANUFACTURED WITHOUT CHLORINE

- Removing chlorine in the mfg. process is more eco-friendly, as this process creates harmful by-products that can be released into the air and water.
- Chlorination can cause gloves to become more brittle over time, reducing their elasticity. It can also leave a residual odor.”



1 SKIN-FRIENDLY ACCELERATOR-FREE

- Manufactured without harmful chemical accelerators - reduces risk of skin irritation.
- Low dermatitis claim recognized by FDA and tested on humans.
- Preventative care to help reduce the risk of developing sensitivity to nitrile gloves (contact dermatitis) due to prolonged glove use.



3 EARTH-FRIENDLY GREEN MANUFACTURING

- Sustainably produced with our innovative, green technology conserving natural resources:
 - 12% less water consumption
 - 17% less CO₂ emissions
- Bright, mint-green color to emphasize the eco-friendly benefits.



4 ENHANCED COMFORT

- Advanced nitrile formulation facilitates optimal ‘touch and feel’ instilling confidence when handling medical instruments.
- Promotes easy donning thanks to an inner polymer coating.



5 CHEMO CLAIM

- Meets USP 800 requirements for PPE compliance for chemotherapy.
- To carry a chemo claim, the glove must adhere to specific requirements set forth by the FDA and pass testing* against permeation of the actual cytotoxic chemicals.

* Testing results available by request.

LET'S SAVE THE WORLD TOGETHER



Our inner boxes are made of 80% recycled material.
Our outer cartons are made of 90% recycled material.



GLOVE the PLANET
www.glovetheplanet.com

Glove the Planet® is our commitment to reducing our environmental impact. By thinking ahead and consuming less, we're able to do more to minimize our global footprint with every step we take. Our commitment to environmental sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations. To learn more, visit www.glovetheplanet.com.



Skin-Friendly Gloves

sempermed® 100+ YEARS
YOUR GLOBAL GLOVE EXPERT™ - EST. 1920 -

Our Gloves are Disposable. The People they Protect are Not.®
Ph 1.800.366.9545 | www.sempermedusa.com
Easy online ordering available on our website **gotgloves**™

REV06092023



LATEX Disposable Gloves



SemperCARE®

FEATURES

- Exam grade latex
- Micro-rough textured surface
- Thickness: Competes with a 5-6 mil glove



*Not applicable where prohibited.

Order # Size

SCLT101	XS
SCLT102	S
SCLT103	M
SCLT104	L
SCLT105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.



StarMed®

FEATURES

- Exam grade latex
- Micro-rough textured surface
- Thickness: Competes with a 4-5 mil glove



*Not applicable where prohibited.

Order # Size

SM101	XS
SM102	S
SM103	M
SM104	L
SM105	XL

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION

84 (14 x 6)

All Sempermed gloves are single use only.

Caution: These products contain natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.

REV05182023



Exam Grade



Aloe Coated



Food Handling Compliant

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920



LATEX Disposable Gloves



SemperGuard® Latex Powder-Free

FEATURES

- Industrial grade latex
- Micro-rough textured surface
- ESD tested and approved
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 4-5 mil glove



*Not applicable where prohibited.

Order # Size

INDPFT101	XS
INDPFT102	S
INDPFT103	M
INDPFT104	L
INDPFT105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
78 (13 x 6)

All Sempermed gloves are single use only.



SemperGuard® Latex Powdered

FEATURES

- Industrial grade latex
- Powdered
- Micro-rough textured surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 4-5 mil glove



*Not applicable where prohibited.

Order # Size

INDPS102	S
INDPS103	M
INDPS104	L
INDPS105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
112 (16 x 7)

All Sempermed gloves are single use only.



GRIPSTRONG® LATEX Powder-Free

FEATURES

- Industrial grade latex
- Micro-rough textured surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove



*Not applicable where prohibited.

Order # Size

GSLF101	XS
GSLF102	S
GSLF103	M
GSLF104	L
GSLF105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
72 (12 x 6)

All Sempermed gloves are single use only.

REV05182023

Caution: These products contain natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.



Industrial Grade



ESD Tested & Approved



Powdered



Food Handling Compliant

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920



VINYL Disposable Gloves



SemperCARE® VINYL STRETCH

FEATURES

- Exam grade vinyl
- Smooth surface
- Stretch properties and cream color
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 4-5 mil glove



Order # Size

SCVNP101	XS
SCVNP102	S
SCVNP103	M
SCVNP104	L
SCVNP105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
84 (14 x 6)

All Sempermed gloves are single use only.



SemperCARE® VINYL Synthetic

FEATURES

- Exam grade vinyl
- Smooth surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove



Order # Size

EVNP101	XS
EVNP102	S
EVNP103	M
EVNP104	L
EVNP105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
84 (14 x 6)

All Sempermed gloves are single use only.



GRIPSTRONG® VINYL

FEATURES

- Industrial grade vinyl
- Smooth surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove



Order # Size

GSVF102	S
GSVF103	M
GSVF104	L
GSVF105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
112 (16 x 7)

All Sempermed gloves are single use only.

REV05182023



Exam Grade



Industrial Grade



Food Handling Compliant

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920



VINYL Disposable Gloves



SemperGuard® Vinyl Powder Free

FEATURES

- Industrial grade vinyl
- Smooth surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove



Order #	Size
VPF102	S
VPF103	M
VPF104	L
VPF105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
96 (16 x 6)

All Sempermed gloves are single use only.



SemperGuard® Vinyl Powdered

FEATURES

- Industrial grade vinyl
- Smooth surface
- Powdered
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3-4 mil glove



Order #	Size
VP102	S
VP103	M
VP104	L
VP105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
112 (16 x 7)

All Sempermed gloves are single use only.



SemperGuard® Black Vinyl

FEATURES

- Industrial grade black vinyl
- Smooth surface
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 3 mil glove



Order #	Size
VBLK102	S
VBLK103	M
VBLK104	L
VBLK105	XL

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

PALLET CONFIGURATION
112 (16 x 7)

All Sempermed gloves are single use only.

REV05182023



Powdered



Food Handling Compliant

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+
YEARS
EST. 1920



POLYETHYLENE Disposable Gloves



SemperGuard® FoodSafe

FEATURES

- Industrial grade polyethylene
- Stretch properties
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 2 mil glove



Order # Size

HPEF202 S

HPEF203 M

HPEF204 L

HPEF205 XL

HPEF206 XXL

GLOVE COUNT

200 Gloves/Box

10 Boxes/Case

PALLET CONFIGURATION

90 (15 x 6)

All Sempermed gloves are single use only.



GRIPSTRONG® POLY

FEATURES

- Industrial grade polyethylene
- Loose fit, easy to don
- Meets FDA CFR 21 Parts 170-199 requirements for food handling
- Thickness: Competes with a 1 mil glove



Order # Size

GSPE504 L

GLOVE COUNT

500 Gloves/Box

20 Boxes/Case

PALLET CONFIGURATION

42 (7 x 6)

All Sempermed gloves are single use only.

REV05182023



Industrial Grade



Food Handling Compliant

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920



SURGICAL Disposable Gloves

sempermed[®] supreme



GLOVE COUNT
50 Pairs/Box ; 6 Boxes/Case

PALLET CONFIG.
32 (8 x 4)

FEATURES

- Sterile Surgical Glove
- Latex
- Micro-rough texture
- Fully anatomical with curved fingers
- Beaded cuff
- AQL: 0.65
- Thickness: Competes with a 4-5 mil glove



Order #	Size
SPFP550	5.5
SPFP600	6
SPFP650	6.5
SPFP700	7
SPFP750	7.5
SPFP800	8
SPFP850	8.5
SPFP900	9

All Sempermed gloves are single use only.



sempermed[®] syntegra cr



GLOVE COUNT
40 Pairs/Box ; 6 Boxes/Case

PALLET CONFIG.
36 (16 x 6)

FEATURES

- Sterile Surgical Glove
- Chloroprene
- Micro-rough texture
- Fully anatomical with curved fingers
- Beaded cuff
- AQL: 0.65
- Thickness: Competes with a 5-6 mil glove

Order #	Size
SCR600	6
SCR650	6.5
SCR700	7
SCR750	7.5
SCR800	8
SCR850	8.5

All Sempermed gloves are single use only.

sempermed[®] syntegra IR



GLOVE COUNT
50 Pairs/Box ; 6 Boxes/Case

PALLET CONFIG.
32 (8 x 4)

FEATURES

- Sterile Surgical Glove
- Polyisoprene
- Micro-rough texture
- Fully anatomical with curved fingers
- Beaded cuff
- AQL: 0.65
- Thickness: Competes with a 4-5 mil glove



Order #	Size
SIR550	5.5
SIR600	6
SIR650	6.5
SIR700	7
SIR750	7.5
SIR800	8
SIR850	8.5
SIR900	9

All Sempermed gloves are single use only.



REV05182023



100+ YEARS

- EST. 1920 -

sempermed[®]
YOUR GLOBAL GLOVE EXPERT[™]

Our Gloves are Disposable.
The People They Protect are Not.[®]

REV05182023

sempermed[®] 
YOUR GLOBAL GLOVE EXPERT[™]

100+ YEARS
- EST. 1920 -

got gloves?

sempermedusa.com | 1-800-366-9545

MILE STONES



1824


GROUND-BREAKING BEGINNINGS

Semperit is founded by Johann Nepomuk Reithoffer, a tailor, entrepreneur and inventor specializing in waterproof fabrics. His patented weaving machine marks the start of the European rubber industry.

OUR FUTURE HOME

Construction commences on a factory in Wimpassing, Austria, which will one day house Sempermed's surgical glove production, research and development center.

1852



1920

FIRST GLOVES

Semperit starts producing latex gloves.

GLOBAL APPEAL

Sempermed surgical gloves are donned by professionals around the world.

1960s



1997

A NEW EUROPEAN MARKET LEADER

Sempermed® Supreme, one of the most frequently used powder-free surgical gloves in Europe, establishes itself as a glove of choice for surgeons.

INNOVATION AWARD

Our patented photochemical process, used to manufacture the skin-friendly, accelerator-free Sempermed® Syntegra UV, is awarded the prestigious EARTO innovation award.

2016



2018

SUSTAINABILITY & SKIN-FRIENDLINESS

The innovative Green Glove Technology introduces a new level of sustainability in glove manufacturing.

A CENTURY OF GLOVE EXPERTISE

2020

100 YEARS

WE DELIVER

Quality & Consistency First®



DEDICATED GLOVE EXPERT

Gloves are our ONLY business. We dipped our first glove 100 years ago. Sempermed has expertise in every phase of the manufacturing process. The result is a product you can trust from people who know their business.

RELIABLE SALES SUPPORT

You can count on your dedicated Account Manager to be there. We are committed to help resolve any issues and find solutions for you and your customers.



DOMESTIC LOGISTICS TEAM



Sempermed's logistics team works with our factories and transportation vendors daily, so you don't have to. SUSA operates in the EST time zone.

CREATIVE MARKETING SUPPORT



Need help with your marketing or branding? Not a problem! Our talented marketing gurus can help with private label artwork and custom marketing pieces tailored to fit your needs. Artwork fees apply.



NO LANGUAGE BARRIER

We communicate with our factories on a daily basis speaking factory "language" for efficient order processing.



IN HOUSE PRODUCT SOURCING

Looking for a specific glove? No problem. Let us find it for you.



DOOR TO DOOR DELIVERY

Saving working capital while a shipment is in transit, terms begin when it lands at your door.

HASSLE FREE FDA & CUSTOMS CLEARANCE



Relax. Let our team of experts get your gloves across the border.

WEEKLY UPDATES

We offer weekly updates on every order you place to keep you informed.



DEDICATED TO QUALITY PRODUCTS, SERVICE & SUPPORT

At Sempermed, we manufacture disposable gloves from start to finish. With more than 100 years of knowledge and experience, we are one of the most efficient and innovative glove manufacturers in the world. We control the process from R&D, production, and logistics to the end user. We are experts at packaging and delivering our full line of products including latex, nitrile, vinyl, chloroprene, and polyethylene gloves. Few companies can promise that kind of dedication to your hand protection. That's why we say Sempermed gloves will meet your expectations, or they will never reach your hands. We are your global glove expert, one that can be trusted unconditionally.



HANDHYGIENE

Infection Prevention

The average
**HEALTH CARE
WORKER HAND**
contains up to
4.6 MILLION
CFUs
(colony forming units)

80%
of illness
causing
germs are
**SPREAD BY
YOUR HANDS**



GLOVES ARE NOT ENOUGH
The World Health Organization (WHO) recommends properly washing hands **BEFORE & AFTER GLOVE USE** by rubbing hands with alcohol-based hand sanitizer or by washing with soap and water



WHEN SHOULD YOU CHANGE YOUR GLOVES?

- 1 As soon as gloves are **DAMAGED**
- 2 **AFTER CONTACT** with a patient or customer and his/her surroundings
- 3 **AFTER CONTACT** with blood, saliva, or other bodily fluids
- 4 When there is evidence of **CONTAMINATION**



PREVENT GERM TRANSMISSION & GLOVE CROSS CONTAMINATION

PROPERLY don & doff gloves	LIMIT the # of surfaces you touch with gloved hands	DISPOSE gloves safely
---	--	------------------------------------




**CHANGE
GLOVES
OFTEN**

On average,
you come in
contact with
**300
SURFACES**
every
30 MINUTES,
exposing
you to
**840,000
GERMS**

HOW OFTEN SHOULD YOU CHANGE YOUR GLOVES?
For maximum protection, it's recommended to change your gloves every:

 CHLOROPRENE 20 - 30 MINS	 NITRILE 30 MINS	 LATEX 15 - 20 MINS	 VINYL 15 MINS OR LESS	 POLYETHYLENE 5-10 MINS OR LESS
---	--	---	---	---

MINIMIZE THE EXPOSURE. MINIMIZE THE RISK.

Let's face it, the world's a pretty germey place. So we've made it our sole mission to provide our customers protection from all those wicked whatever's. Since we focus only on manufacturing and supplying gloves, we set the bar on quality standards that surpass regulatory guidelines for both consistency and performance. Because there's a person attached to every single one of the hands we protect - we feel that it's people who we're really protecting. We take enormous pride in being the first and last line of defense for people in harm's way.

Source: World Health Organization, 2009.

sempermed[®]
YOUR GLOBAL GLOVE EXPERT™

100+
YEARS
EST. 1920


Our Gloves are Disposable.
The People they Protect are Not.[®]
www.sempermedusa.com • Ph 1.800.366.9545

Rev10012021



FACING FENTANYL

Stay Informed. Stay Protected.



FENTANYL WAS INVENTED IN 1960
for treating severe pain due to illnesses such as cancer or as a surgical anesthetic




Fentanyl is often **ACCIDENTALLY ABUSED** & commonly mixed with heroin or cocaine in **POWDER, PILL, OR LIQUID FORM** and is often manufactured in China.



Fentanyl is a **HIGHLY ADDICTIVE** synthetic opioid that is **80 TO 100 TIMES** more potent than morphine

FENTANYL IS SO POTENT that a dose the size of a **GRAIN OF SAND** can be deadly



YOU CAN'T SEE, SMELL OR TASTE FENTANYL




Fentanyl can easily be **ABSORBED THROUGH THE SKIN**, putting medical staff, first responders, police, security, coroners, and morticians **AT RISK**



For the first time in US history, our nation surpassed **100,000 OVERDOSE FATALITIES** in a 12 month period.

Fentanyl has quickly & quietly **BECOME THE LEADING CAUSE OF DEATH** among adults **AGES 18-45.**



80% of Fentanyl in America comes into our country **THROUGH THE US-MEXICO BORDER**



OVER 64,000 of those deaths **WERE DUE TO FENTANYL.**



The CDC confirms that Fentanyl is now the **DEADLIEST DRUG IN AMERICA**

For more information about protection against fentanyl, please contact your sales rep.



got FENTANYL gloves?

As a global glove expert, it is our job to educate about the risks you may face on the job, such as Fentanyl exposure. To help minimize the risk, Sempermed offers a handful nitrile exam gloves that provide protection against harmful fentanyl contact. Choose from SemperShield® Standard Cuff, SemperShield® Extended Cuff, SemperForce® Orange, SemperForce® Blue, SemperForce® Black, SemperCare® Nitrile, StarMed® PLUS, StarMed® Ultra, StarMed® Select, Tender Touch® 100, Tender Touch® 200, SemperSure® and SemperSure® Black.

Testing for protection against Fentanyl is conducted by an independent lab and verified by SUSA Regulatory. Testing results are available by request.

Sources: <http://www.murphy.house.gov>, <http://www.fda.gov/oc/fentanyl-overdose>



Our Gloves are Disposable. The People they Protect are Not.
www.sempermedusa.com • Ph 1.800.366.9545




CHEMOGLOVES

Stay Informed. Stay Protected.




According to the World Health Organization, there were **14.1 MILLION NEW CANCER CASES & 32.8 MILLION PEOPLE LIVING WITH CANCER** (within 5 years of diagnosis) in 2012 worldwide


At present, more than **100** ANTICANCER DRUGS are used in daily clinical practice



Our Chemotherapy tested gloves meet **USP 800** REQUIREMENTS FOR PPE COMPLIANCE FOR CHEMOTHERAPY



In order for a glove to carry a chemo claim, the glove must adhere to specific requirements set forth by the **FDA** & PASS TESTING AGAINST PERMEATION of the actual cytotoxic chemicals



The **EXCRETIONS** OF PATIENTS such as urine, feces, vomit, and sweat, may be **CONTAMINATED** with chemo drugs, making urine bags, clothing and hospital beds hazardous objects



OUR GLOVES ARE DISPOSABLE. THE PEOPLE THEY PROTECT ARE NOT®

We have made it our mission to provide our customers with protection from all those wicked whatever's. Since there is a person attached to every single one of the hands we protect, we feel it is the people themselves that we are really protecting. If you handle cytotoxic chemicals or work in oncology, try one of our nitrile exam glove brands that have been tested for use against chemotherapy drugs*: SemperShield® Extended Cuff, SemperSure® Black, SemperSure® Green, SemperSure® Nitrile, SemperForce Blue®, SemperCare®, TenderTouch®, TenderTouch® 100, StarMed® Plus, StarMed® Ultra, and StarMed® Select.

*Meets USP 800 requirements for PPE compliance for chemotherapy



ESD GLOVES

We've Got You Covered



WHAT IS ESD

ELECTROSTATIC DISCHARGE:

The sudden flow of electricity between two electrically charged objects caused by contact, an electrical short, or dielectric breakdown

STATIC ELECTRICITY CAN POSE A FIRE-RISK



In industries that handle flammable liquids, such as automotive, painting, chemical, lab, pyrotechnic and petroleum/gas



The human body is an **EXCELLENT CONDUCTOR OF ELECTRICITY** and a person can become statically charged with upwards of **SEVERAL THOUSAND VOLTS**




The static charge that builds up in humans must be released, usually resulting in **AN ELECTRIC SHOCK**

By definition, a static dissipative material has a **SURFACE RESISTIVITY** of 1×10^6 to 1×10^{12} ohms per square

ESD GLOVES

provide **HAND PROTECTION** that **HELPS PREVENT PRODUCT DAMAGE** from electrostatic discharge and are **GENERALLY SAFE** to work with flammable liquids/chemicals or static sensitive devices



In many industries, this discharge of static electricity can **DAMAGE A VARIETY OF PRODUCTS**, especially sensitive electronics



ESD TESTED & APPROVED GLOVES

Are you handling sensitive electronics or flammable liquids? We've got you covered! Sempermed offers ten brands that are ESD tested and approved. ESD gloves provide hand protection that helps prevent product damage from electrostatic discharge, and they are safe to work with static sensitive devices where bare hands are not desired*. We can send you free samples of our ESD tested and approved gloves: SemperSure® Black Nitrile, SemperForce® Black Nitrile, SemperForce® Blue Nitrile, SemperCare® Nitrile, SemperShield® standard cuff, SemperGuard® Nitrile, SemperGuard® HD Nitrile, GripStrong® Nitrile, GripStrong® White Nitrile, and SemperGuard® Latex.

* Level of ESD protection should be confirmed for your needs.



SENSITIVE SKIN?

We've Got You Covered




SENSITIVITY TO LATEX PROTEINS
remains a primary cause of concern among glove users

LATEX PROTEIN ALLERGIES AFFECT
4.3%
of the general population




ALOE
has a soothing effect on the skin and has been shown to **HELP REDUCE ITCHINESS & INFLAMMATION**

TYPE I (IMMEDIATE)
Hypersensitivity to natural rubber latex proteins can be **VERY SERIOUS OR OCCASIONALLY EVEN FATAL**



TYPE IV (DELAYED)
Hypersensitivity which involves a skin rash and is often **CAUSED BY CHEMICALS** used in the manufacturing process




SOME GLOVE USERS can develop rashes linked to irritant contact **DERMATITIS**. This may be caused by chemicals such as **ACCELERATORS** used in the manufacture of the gloves




Many hospitals and healthcare facilities have **MIGRATED AWAY FROM LATEX GLOVES**, yet some users experience skin sensitivities to non-latex gloves

TRY ACCELERATOR-FREE GLOVES TO MINIMIZE THE RISK



ARE YOU SENSITIVE TO YOUR GLOVES?

Do you have delicate skin or suffer from dry, chapped hands? Do you often wonder if you are sensitive to your gloves? We've got you covered! Sempermed offers a handful of glove brands that provide relief for irritated hands. Best Touch® StarMed® ROSE nitrile exam has an aloe coating, which soothe and nourish skin irritation. Our SemperSure® glove family are accelerator-free, low-dermatitis, nitrile exam gloves designed for those with skin sensitivities to standard nitrile gloves. They offer all the benefits of a nitrile glove, minus the chemical accelerators (added in the manufacturing process) that may cause skin reactions.



PROP 65

Understanding DOTP & DINP Gloves

ENACTED IN 1986,
PROP 65
IS A CALIFORNIA LAW



intended to reduce or eliminate exposure to certain substances. It is comprised of a list of substances that the State of California recognizes as causing cancer and/or birth defects

DINP

(Diisononyl phthalate) is a plasticizer used in the manufacture of **MANY VINYL** (polyvinyl chloride/PVC) **PRODUCTS**, including electrical wire, coated fabrics, automotive parts, vinyl flooring, & footwear, as well as general purpose & medical examination **GLOVES**



in 2013 the State of California
ADDED
DINP
to the
Proposition 65
list of chemicals
"known to the
State of California
to cause cancer"



The addition of DINP to the Prop 65 list of banned chemicals is forcing glove manufacturers to re-formulate to

**A NON-BANNED
PLASTICIZER CALLED**

DOTP

PROP 65 LABEL WARNING

As of August 2018

Clear & reasonable label warning requirements include approved "short form" & "long form" versions



WARNING: Cancer -
www.P65Warnings.ca.gov

ALL
Sempermed®

VINYL GLOVES ARE

100%

FREE

OF DEHP & DINP

SEMPERMED DOTP VINYL GLOVES

There is scientific and medical evidence that vinyl gloves are safe and serve a valuable role in the medical and industrial glove market. Simply because DEHP & DINP have been placed on the Prop 65 list, does NOT mean that vinyl gloves containing these substances are unsafe or pose health risks. Nonetheless, all Sempermed USA vinyl gloves are DOTP to avoid mandatory Prop 65 warning labeling requirements. Try free samples of any of our vinyl gloves: SemperCare® Stretch Vinyl, SemperCare® Synthetic Vinyl, GripStrong® Vinyl, SemperGuard® (Powdered and Powder-Free) Vinyl, and SemperGuard® Black Vinyl.



gotgloves

PRODUCT
CATALOG

Our Gloves are Disposable.
The People they Protect are Not.®

www.sempermedusa.com • Ph 1.800.366.9545



Rev 12/20/2022

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+
YEARS
EST. 1920

UDI LABELING

Enforcement & Data Submission

September 24, 2013 the



**ESTABLISHED
A UNIQUE
IDENTIFICATION
SYSTEM**
to adequately
identify medical devices

GTIN
GLOBAL TRADE ITEM
NUMBER

an item barcoding system
developed by the
non-profit organization
GS1. GS1 IS AN
FDA-ACCREDITED
ISSUING AGENCY

UDI

UNIQUE DEVICE
IDENTIFICATION #

Identifies a product. Always
numeric, non-intelligent
number, unambiguous



ENFORCEMENT OF UDI
LABELING and data submission
for Class 1 medical devices,
**INCLUDING DISPOSABLE
EXAMINATION AND
SURGICAL GLOVES,**
has been extended again to



To be compliant, the labeler must
submit information of each device to
FDA's Global Unique Device
Identifier Database (GUDID), which
serves as the repository of key
device identification data. **THE
LABELER IS DEFINED AS THE
OWNER OF THE BRAND. THE
OWNER OF A PRIVATE LABEL
BRAND IS THEREFORE THE
LABELER.** All such devices must be
properly labeled with a UDI and be
registered on GUDID by
DECEMBER 8, 2022
by the labeler

GS1 develops global
standards for the
IDENTIFICATION
OF GOODS



UPC
UNIVERSAL
PRODUCT CODE

NEW GTIN UPC
(inner boxes)



12 DIGITS

NEW GTIN UPC
(outer cartons)



14 DIGITS

ENFORCEMENT OF UDI LABELING

The FDA has again extended the deadline for the enforcement of UDI requirements for Class 1 medical devices, including disposable examination and surgical gloves, to 8 December 2022. To meet this compliance date, Sempermed has labeled each of our products with a UDI, in the form of 12 and 14 digit Universal Product Codes (UPC), which were created using GTIN (barcoding) software provided by GS1, and registered on GUDID. Accurate and synchronized data improves supply chain efficiency and patient safety by reducing picking errors, facilitating faster product recalls, enabling efficient tracking of medical devices, while enhancing inventory management.



GIVING BACK

Sempermed Proudly Supports



Clothes To Kids



What was the inspiration for the mission behind Clothes To Kids?

\$65 = 1 WARDROBE

Since 2003, Sempermed has donated more than

4,054

WARDROBES

A WARDROBE INCLUDES:

- 5 pairs of new underwear
- 5 pairs of new socks
- 5 tops (including uniforms)
- 4 bottoms (including uniforms)
- 1 dress (optional)
- 1 pair of shoes
- 1 jacket (in season)



THE HIGHEST RATING given to a charity

IT ALL STARTED WITH ONE PURPLE COAT
Read the full story by visiting www.clothestokids.org

Clothes To Kids envisions a community in which EVERY SCHOOL-AGE CHILD HAS QUALITY CLOTHING so that they may attend school with the CONFIDENCE & SELF-ESTEEM needed to achieve academic success



89%

OF EXPENDITURES go directly to the program

2023 ANNUAL CHARITY EVENTS



APRIL 20
6PM-9PM

21st Birthday Celebration
The Coliseum - St. Petersburg, FL
An evening of delicious food, open bar, live music, dancing, and fundraising



JULY & AUGUST
(SUSA COLLECTS IN MAY)

Drop Off Your Drawers
New Underwear Drive



OCTOBER

Shoetober - New Shoe Drive & Kick One Back Shoe-raiser Event

GRASS ROOTS

3 store locations in Florida
& affiliate locations in
CO, NC, CT, & RI

CLOTHE A KID. CHANGE A LIFE.

Clothes To Kids was founded in 2002 in Pinellas County, Florida, home of Sempermed USA's headquarters. Their mission is to provide new & used quality clothing to low-income or in crisis school-age children, free of charge. \$65 will provide one child with a wardrobe; that's a full week's worth of school clothes. **Want to help? Visit www.clothestokids.org. All donations are tax-deductable.**

sempermed[®]
YOUR GLOBAL GLOVE EXPERT™

100+
YEARS
EST. 1922

Our Gloves are Disposable.
The People they Protect are Not.[®]
www.sempermedusa.com • Ph 1.800.366.9545



Rev06202023

GIVING BACK

Sempermed Supports



Prostate Cancer
Foundation
Curing Together.



ABOUT 1 IN 8 MEN WILL BE DIAGNOSED with Prostate Cancer during his lifetime



Prostate Cancer is the **MOST COMMON** type of cancer found in American men (after skin cancer)

A MAN GETS DIAGNOSED with Prostate Cancer every



ALL MEN should consider starting routine screening for Prostate Cancer around the age **OF 50**

The PCF has helped reduce the death rate of Prostate Cancer in the United States by

52%



99% of patients diagnosed with local or regional disease live **5 YEARS** OR LONGER AFTER DIAGNOSIS

Black men are about **75% MORE LIKELY** to develop Prostate Cancer than non-Hispanic White men, and are more than twice as likely to die from the disease.



The Prostate Cancer Foundation has raised **NEARLY \$1 BILLION** for prostate cancer research (since 1993) to **SAVE & EXTEND MEN'S LIVES**



EVERY CASE LENDS A HAND AGAINST PROSTATE CANCER

Sempermed is committed to supporting the Prostate Cancer Foundation (PCF) and its efforts in research and awareness of the disease. As an organization, it is important for us to donate to a cause where funds are directly applied towards developing a cure. That's why Sempermed is donating a portion of the proceeds from each case of StarMed® PLUS, SemperForce®Black, and Best Touch® StarMed® ROSE nitrile exam gloves to the Prostate Cancer Foundation. Scan the QR codes below to learn more about each product.



SMNP PRODUCT INFO



BKMF PRODUCT INFO



SMNR PRODUCT INFO

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920

Rev08152022

SUSTAINABILITY

Contributing to a Prosperous world

COMPLIANCE & INTEGRITY

are the foundation of our sustainability strategy
SEMPERMED'S
5 CORE PILLARS
OF SUSTAINABILITY:

1. RESOURCE MGMT AND ENVIRONMENT
2. HEALTH AND SAFETY
3. EMPLOYEES AND SOCIETY
4. SUPPLIERS
5. INNOVATION



SAFETY FIRST

We strive to PREVENT ACCIDENTS and continuously IMPROVE WORKING CONDITIONS to retain our highly qualified, enthusiastic workforce. To achieve this, we remain focused on **TRANSPARENCY, INTEGRITY AND CONSISTENCY**



Business Social Compliance Initiative

The Semperit Group has committed itself to high ethical and social standards, as well as fair working conditions for all employees world-wide. We are a proud BSCI member and proponent of employee rights.



HEALTH, SAFETY, ENVIRONMENT, ENERGY and QUALITY (HSEQ) are integral parts of our business.



We focus on encouraging all employees to achieve continuous improvement and are committed to our long-term goal of **ZERO INCIDENTS**. We strive to **AVOID WASTE** and higher consumption in order to ensure **RESOURCE-EFFICIENT, SUSTAINABLE MANUFACTURING OF FIRST-CLASS PRODUCTS**.

The objective of the EcoVadis Corporate Social Responsibility (CSR) rating is to measure the quality of a company's CSR management system through its policies, actions, and results.



2022: Semperit Group has achieved the **PLATINUM STATUS ECOVADIS CSR RATING** - Top 1% of all companies rated



2021: Semperit Group has achieved the **GOLD STATUS ECOVADIS CSR RATING**.



2020 & 2019: Our Wimpassing site was awarded the **GOLD STATUS ECOVADIS CSR RATING** for performance in the area of corporate social responsibility.

THIS SEAL IS OUR PROMISE FOR:

- Sustainably produced skin-friendly gloves
- Strict quality controls for optimal protection
- Reliability by ensuring European standards
- Use of renewable energies in the production process
- Innovative mfg. processes and products



YOU CAN'T MANAGE WHAT YOU DON'T MEASURE

Semperit participates in the international CARBON DISCLOSURE PROJECT (CDP) emissions reporting program, joining over 9,600 other companies who report on CLIMATE CHANGE, SUPPLY CHAIN, WATER SECURITY & FORESTS.

COMMITTED TO SUSTAINABLE GROWTH

We are aware of our responsibility towards customers, employees, suppliers, shareholders, society, and the environment. We handle our resources carefully, treat our fellow human beings with respect, and consider the consequences of our actions at all times. That is why we have anchored sustainability within our strategic corporate goals and formulated correspondingly strict standards. Our group-wide sustainability policy is based on our code of conduct. This is supplemented by globally applicable guidelines and standards in the five focus areas - resource management, health and safety, work and society, suppliers, and innovation. We are taking clear steps towards a sustainable company development. We can make a difference together.

That's how we GLOVE the PLANET®

GLOVE THE PLANET® IS OUR COMMITMENT TO REDUCING OUR ENVIRONMENTAL IMPACT.

By thinking ahead and consuming less, we're able to do more to minimize our global footprint with every step we take. Our commitment to environmental sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations.



Recycled Packaging

We Save Trees!

- We use over **5 million pounds of recycled materials** in our packaging every year.
- Our recycled packaging **saves approximately 5,000 trees** annually.



Reduce Waste

Taking Steps in the Right Direction

- Our goal is to **reduce our material waste by 4%** each year.
- We produce flip flops from **leftover nitrile glove scraps**.



Green Manufacturing

We Save Resources

- Using no chlorine and **13% less water and gas** in our innovative green mfg. process, reduces our **CO₂ impact by more than 17%**.



30 by 2030

We Save Energy

- We have a mission to **reduce** our overall energy, emissions, waste, and water **by 30%** by the year 2030. That's a **3% reduction each year**.



got gloves?



GLOVE the PLANET

SCAN ME



SCAN TO LEARN MORE ABOUT OUR COMMITMENT TO SUSTAINABILITY

Rev01242022

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920

Our Gloves are Disposable. The People they Protect are Not.®

www.sempermedusa.com • Ph 1.800.366.9545





GLOVE *the* **PLANET**[®]

www.GlovethePlanet.com

OUR COMMITMENT TO REDUCING OUR ENVIRONMENTAL IMPACT

By thinking ahead and consuming less, we're able to do more to minimize our global footprint with every step we take. Our commitment to environmental sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations.

RECYCLABLE

Alternative to conventional waste disposal by converting uncontaminated waste materials into new materials.

Non-contaminated recycled gloves can be up-cycled into new eco-friendly plastic products.



COMPOSTABLE

Compostable plastics are those that biodegrade in industrial composting operations at the rate of other compostable materials.

Compostable gloves eventually break down into dirt.



BIODEGRADABLE

All organic matter will eventually biodegrade. The term biodegradable must be qualified by the environment and timeframe.

Biodegradable gloves eventually break down into small pieces.

Biodegradable does not mean a material is compostable or recyclable.



Only Non-contaminated gloves can be recycled or sent to the land fill. Gloves in contact with blood, body fluids or chemicals must be disposed of as waste.



LATEX GLOVES

Can take 1-5 Years* to break down



NITRILE GLOVES

Can take Decades* to break down



VINYL GLOVES

Can take over 100 years* to break down



BIODEGRADABLE GLOVES

Can take up to 5 years* to break down

*Depends on landfill conditions such as humidity & temperature



RAW MATERIALS ARE LIMITED

Since 1970 the global raw material consumption has more than tripled.

ENERGY IS VALUABLE

By 2050 the global energy demand will increase by 70%-100%.

CLIMATE CHANGE

Is thought to be the #1 threat to global public health today.

Sources: Definitions: www.merriam-webster.com, www.businessdictionary.com Department of Ecology (WA), Focus on "Biobased," "Biodegradable," & "Compostable" Plastics , Aug. 2014



SCAN ME




Our Gloves are Disposable. The People they Protect are Not.[®]
Ph 1.800.366.9545 | www.sempermedusa.com
www.glovetheplanet.com



REV02062023

got gloves?



QUALITY GLOVES BY **sempermed**[®] THAT LEAVE A **LASTING IMPRESSION**

ENHANCED COMFORT & FIT • SUPERIOR PROTECTION • LESS HAND FATIGUE



SEMPERFORCE **BLACK** Nitrile Exam

REORDER# **BKNF102-106**

Sizes **S - XXL**

- Strong & durable
- Excellent gripping capabilities

100 gloves/bx (90 XXL) | 10 bxs/cs



Best Touch StarMed **ROSE** Nitrile Exam with Aloe

REORDER# **SMNR201-205**

SIZES **S - XL**

- Coated w/ Aloe to hydrate dry hands
- Helps keep cuticles soft

200 gloves/bx (180 XL) | 10 bxs/cs



SemperSure **Black** Nitrile Exam

REORDER# **SUNB102-105**

SIZES **S - XL**

- Accelerator-free/low-derm
- Ideal for sensitive skin

100 gloves/bx | 10 bxs/cs



LOVE the PLANET
www.glovetheplanet.com

Glove the Planet[®] is our commitment to reducing our environmental impact. By thinking ahead and consuming less, we're able to do more to minimize our global footprint with every step we take. Our commitment to environmental sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations. To learn more, visit www.glovetheplanet.com.

PRODUCT CATALOG



sempermed[®]
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
EST. 1920

Our Gloves are Disposable. The People they Protect are Not.
Ph 1.800.366.9545 | www.sempermedusa.com
Easy online ordering available on our website **gotgloves**SM

REV05182023