

black nitrile examination gloves

S^{semper} sapphire[™] XP black

Semperguard[®] sapphire[™] XP black exam nitrile gloves are designed to offer excellent grip and tactile sensitivity, making them ideal for a variety of tasks that require precision and control. The black color of the gloves are especially appealing in industries such as beauty, laboratory, food services, janitorial/ sanitation, and safety. The dark color helps conceal color transfer, stains, and other messes, maintaining a cleaner and more professional appearance throughout use. This feature is especially useful in environments where aesthetics and cleanliness are a priority.



- Sleek and professional appearance Hides stains caused from grease, dirt, dyes, and ink
- Textured fingertips
 Enhanced dexterity and gripping capabilities
- Tested against 33 chemotherapy drugs, fentanyl, and gastric acid Meets USP 800 requirements for PPE compliance for chemotherapy
- Safe for food contact applications Meets FDA 21 CFR Part 177 requirements for food handling
- Made in Malaysia
 Manufactured in Malaysia using stateof-the-art technology to achieve quality that exceeds FDA standards

100 Gloves/Box | 10 Boxes/Case (90 XXL)



Semper Sapphire[™] guard XPblack

black nitrile examination gloves

| Technical Data | | | | | |
|--------------------------------------|--|--|--|--|--|
| ТҮРЕ | Nitrile single use disposable protective glove, non-sterile | | | | |
| GRADE | Exam | | | | |
| COLOR | Black | | | | |
| MATERIAL | Nitrile Butadiene Rubber (NBR) not made with natural rubber latex | | | | |
| GLOVE SURFACE | Finger textured | | | | |
| COUNTRY OF ORIGIN | Malaysia | | | | |
| SHELF LIFE | 3 years in original packaging if stored as per DIN 7716, ISO 2230 | | | | |
| SIZE | Length: 240mm (min) Width: 95mm +/- 10mm (based on size M) | | | | |
| WALL THICKNESS | Finger Palm .11mm .07mm 4.3mils 2.8mils based on average production value measured in single layer | | | | |
| TENSILE STRENGTH (after aging) | 25MPa (min) | | | | |
| ELONGATION (after aging) | 400% (min) | | | | |
| BARRIER PERFORMANCE | AQL ≤ 1.0 | | | | |
| QUALITY CERTIFICATION | amfori-BSCI, RGA, Sedex, WRAP | | | | |
| TESTING | ASTM D412, ASTM D573, ASTM D3767, ASTM D5151, ASTM D6124, ASTM D6978, ASTM F1671/1671M | | | | |

For maximum protection, we recommend changing nitrile gloves every 20-30 minutes



PEOPLE, PROTECTING PEOPLE

Packaging Specifications

DIMENSIONS

Dispenser 245L x 125W x 42H mm Carton 260L x 223W x 255H mm

Recommended Applications



| Reorder Numbers / Glove Counts* | | | | | | |
|---------------------------------|----------|----------|---------------|-----------|--|--|
| | GSBNF101 | X-Small | 100 gloves/bx | 10 bxs/cs | | |
| | GSBNF102 | Small | 100 gloves/bx | 10 bxs/cs | | |
| | GSBNF103 | Medium | 100 gloves/bx | 10 bxs/cs | | |
| | GSBNF104 | Large | 100 gloves/bx | 10 bxs/cs | | |
| | GSBNF105 | X-Large | 100 gloves/bx | 10 bxs/cs | | |
| | GSBNF106 | XX-Large | 90 gloves/bx | 10 bxs/cs | | |

SAMPLES: Available by request - ask your sales rep! *glove counts "by weight"



04222025USA

INPORTANT NOTE: The latest product information is available at www.gotgloves.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. HARPS is not liable for incorrect use of the gloves nor will HARPS be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used. In case of advice should be sought before use. The information and leasification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. HARPS bears no liability for allergic reactions. @Copyright 2024 HARPS



www.gotgloves.com | 800.366.9545