





# coated with aloe



Sempermed® rose™ exam nitrile gloves are designed to elevate comfort and protection in any setting, from exam rooms to tattoo shops to day spa. These gloves feature a soothing aloe coating that helps prevent skin irritation, keeping hands moisturized and comfortable during extended use. The unique formulation of these gloves enhances protection while maintaining exceptional tactile sensitivity, allowing for precise handling without compromising comfort.

- Coated with aloe: Helps moisturize dry hands and helps prevent skin irritation while you work. A great choice for sensitive skin.
- Fun, unique color: Unique fuchsia pigment color brightens any environment from the ER to the day spa.
- Protects users and patients from the transmission of infections: Designed to provide protection and barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.

#### **Technical data**

TYPE	Nitrile single use disposable protective glove, non-sterile		
GRADE	Examination		
COLOR	Fuchsia		
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex		
GLOVE SURFACE	Finger textured		
COUNTRY OF ORIGIN	Malaysia		
SHELF LIFE	3 years in original packaging if stored as per DIN 7716, ISO 2230		
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)		
WALL THICKNESS	Finger Palm .07mm .05mm 2.8mils 2.0mils based on average production value measured in single layer		
TENSILE STRENGTH (after aging)	14MPa (min)		
ELONGATION (after aging)	400% (min)		
BARRIER PERFORMANCE	AQL <u>&lt;</u> 1.5		
TESTING	ASTM D412, ASTM D573, ASTM D3767, ASTM D5151 and ASTM D6124		

### **Purpose**

Single use medical device class 1, as per the FDA Single-use medical device

Suitable for food handling, as per relevant sections of 21 CFR Parts 170-199

## **Packaging**

DIMENSIONS

Dispenser 245L x 125W x 74H mm Carton 385L x 255W x 255H mm

#### Reorder numbers / Glove counts\*

	SMNR201	X-Small	200 gloves/bx	10 bxs/cs
	SMNR202	Small	200 gloves/bx	10 bxs/cs
	SMNR203	Medium	200 gloves/bx	10 bxs/cs
Ī	SMNR204	Large	200 gloves/bx	10 bxs/cs
	SMNR205	X-Large	180 gloves/bx	10 bxs/cs

SAMPLES: Available by request - ask your sales rep!

For maximum protection, we recommend changing nitrile gloves every 20-30 minutes "glove counts "by weight"

