







sempermed® velvet™ is a lightweight and soft examination and protective glove. It is a facility-wide glove and is suitable for a wide range of applications. This glove is the perfect latex-free alternative for sensitive skin. The soft nitrile formulation allows the glove to conform to the hand like a second skin, combining optimal fit and comfort with unrestricted hand motion.

- Soft as velvet: the soft formulation of the glove offers enhanced comfort and fit.
- Your daily choice: suitable for a wide variety of medical applications where protection is critical.
- Tested for protection against 56 chemotherapy drugs, fentanyl, and gastric acid.
- Protects users and patients from the transmission of infections. Designed to provide a high level of protection and barrier against contaminants, bacteria, and viruses, reducing the risk of crosscontamination and infection in a clinical setting.

Technical data

TYPE	Nitrile single use disposable protective glove, non-sterile		
GRADE	Examination		
COLOR	Violet-Blue		
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex		
GLOVE SURFACE	Finger textured		
COUNTRY OF ORIGIN	Malaysia		
SHELFLIFE	3 years in original packaging if stored as per DIN 7716, ISO 2230		
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)		
WALL THICKNESS	Finger Palm .09mm .06mm 3.6mils 2.4mils based on average production value measured in single layer		
TENSILE STRENGTH (after aging)	28MPa (min)		
ELONGATION (after aging)	400% (min)		
BARRIER PERFORMANCE	AQL <u><</u> 1.0		
QUALITY CERTIFICATION	amfori-BSCI, RGA, Sedex, WRAP		
TESTING	ASTM D412, ASTM D573, ASTM D3767, ASTM D5151, ASTM D6124, ASTM D6978, and ASTM F1671/1671M		

Purpose

Single use medical device class 1, as per the FDA Single-use medical device

Single use medical examination glove, as per ISO 11193

Suitable for food handling, as per relevant sections of 21 CFR Parts 170-199

Packaging

300ct	Dispenser 245L x 128W x 83H mm
DIMENSIONS	Carton 430L x 266W x 257H mm
100ct	Dispenser 245L x 125W x 42H mm
DIMENSIONS	Carton 260L x 223W x 255H mm

Reorder numbers / Glove Counts*

SMNP301		300 gloves/bx	
SMNS101	X-Small	100 gloves/bx	10 bxs/cs
SMNP302 SMNS102	Small	300 gloves/bx 100 gloves/bx	10 bxs/cs
SMNP303 SMNS103	Medium	300 gloves/bx 100 gloves/bx	10 bxs/cs
SMNP304 SMNS104	Large	300 gloves/bx 100 gloves/bx	10 bxs/cs
SMNP305 SMNS105	X-Large	250 gloves/bx 100 gloves/bx	10 bxs/cs

SAMPLES: Available by request - ask your sales rep! *glove counts "by weight" For maximum protection, we recommend changing nitrile gloves every 20-30 minutes









ev 04/2025 US